


Model 96A
V/f Preamplifier

Part No. 617282 Rev. B

Copyright © 1997
Physical Electronics, Inc.
6509 Flying Cloud Drive
Eden Prairie, MN 55344

The PHI logo is a registered trademark of Physical Electronics, Inc.

Physical Electronics, PHI, and  are trademarks of Physical Electronics, Inc.

Limited Warranty

Except as otherwise provided herein, the Seller warrants to Buyer that the equipment sold hereunder, whether it is new equipment or remanufactured (reconditioned) equipment, is, at the time of shipment to Buyer from Seller, free from defects in material and workmanship. As Buyer's sole exclusive remedy under this warranty Seller agrees either to repair or replace, at Seller's sole option and free of part charge to Buyer, any part or parts of such equipment which, under proper and normal conditions of use prove to be defective within 12 months from the date of receipt by the Buyer. Warranty period for equipment requiring installation by Seller will commence on completion of standard installation services. If customer delays installation beyond 45 days after delivery, the warranty period will commence to run 45 days after delivery. After installation, any realignment, readjustment, recleaning or recalibration, provided it does not relate to a proven defect in material or workmanship, shall be performed only at Seller's then current rates for service.

Exclusions and Limitations

It is recognized that some parts by their nature (expendable items) may not function for one full year; therefore, excluded from the foregoing warranty are filaments, anodes, cathodes, multipliers, retard grids, special ceramics, ionizers, along with other such parts mentioned in the applicable operating manual.

The foregoing warranty excludes certain major items or accessories specifically indicated on applicable price lists or quotations, as to which Seller passes to Buyer whatever warranty is provided to Seller by the manufacturer or the

specific warranty indicated by the price list or quotation.

This warranty does not cover loss, damage, or defects resulting from transportation to the Buyer's facility, improper or inadequate maintenance by Buyer, buyer-supplied software or interfacing, unauthorized modification or misuse, operation outside of the environmental specifications for the equipment or improper site preparation and maintenance.

Product Service

All claims must be brought to the attention of Seller within 30 days of the failure to perform.

Seller at his option may require the product to be returned to the factory, transportation prepaid, for repair.

Refund of Purchase Price

In lieu of the foregoing, Seller may at any time elect, in its sole discretion, to discharge its warranty by accepting the return of such equipment and refunding any portion of the purchase price paid by Buyer.

Software and Firmware Products

The sole exclusive warranty applicable to software and firmware products provided by Seller for use with a processor will be as follows: Seller warrants that such software and firmware will conform to Seller's program manuals current at the time of shipment to Buyer when properly installed on that processor. Seller does not warrant that the operation of the processor software or firmware will be uninterrupted or error free.

No other warranty is expressed or implied. Seller expressly disclaims the implied warranties of merchantability and fitness for a particular purpose.

OPERATOR SAFETY SUMMARY

All PHI Systems have been designed to assure operator safety. However, like all other sophisticated instruments, continual operator safety is dependent on the proper use of system components. Such use is described in the manuals supplied with each unit.

LISTED BELOW ARE WARNINGS APPLICABLE TO THIS EQUIPMENT. ALL PERSONNEL INVOLVED IN THE OPERATION AND MAINTENANCE OF THIS EQUIPMENT MUST FULLY UNDERSTAND THE WARNINGS AND THE PROCEDURES BY WHICH THE HAZARD IS TO BE REDUCED OR ELIMINATED.

WARNING

THE PRODUCT(S) COVERED IN THIS MANUAL HAS BEEN DESIGNED TO FUNCTION SAFELY WITH THE REQUIRED DEGREE OF PRECISION WHEN USED IN THE PRESCRIBED MANNER.

WE DO NOT RECOMMEND THAT THIS EQUIPMENT BE MODIFIED FOR ANY NON-STANDARD APPLICATION SINCE HAZARDOUS CONDITIONS MAY RESULT. DUE TO THE FACT THAT PHYSICAL ELECTRONICS DIVISION OF PERKIN-ELMER HAS NO CONTROL OVER CUSTOMER MODIFICATIONS TO PHI PRODUCTS SHIPPED, IT DISCLAIMS ANY RESPONSIBILITY FOR ANY MALFUNCTIONS OR ACCIDENTS THAT MAY RESULT!

DANGER

ELECTRICAL SHOCK HAZARD

HIGH VOLTAGES ARE PRESENT IN THE SYSTEM WHEN THE SYSTEM POWER INPUT LINES ARE CONNECTED. DISCONNECT INPUT POWER AT THE WALL BEFORE MAKING ANY ADJUSTMENTS. REFER SERVICING TO PERSONNEL WHO HAVE BEEN TRAINED AND HAVE WORKING EXPERIENCE WITH VOLTAGES IN EXCESS OF 50 VOLTS.

ALL ELECTRICAL CABLES ASSOCIATED WITH VARIOUS UNITS INCLUDED IN A SYSTEM ARE WELL SHIELDED. HOWEVER, CARE MUST BE TAKEN NEVER TO COME IN CONTACT WITH ANY ASSOCIATED TERMINALS WHEN THE POWER IS ON. SOME OF THESE LEADS CARRY POTENTIALLY LETHAL HIGH VOLTAGES. OTHER LEADS MAY CARRY SUFFICIENT RF POWER TO INFLICT SEVERE BURNS.

RF INTERFERENCE

THIS EQUIPMENT GENERATES, USES, AND CAN RADIATE RADIO-FREQUENCY ENERGY, AND IF NOT INSTALLED AND USED IN ACCORDANCE WITH THE INSTRUCTION MANUAL, MAY CAUSE INTERFERENCE TO RADIO COMMUNICATIONS. OPERATION OF THIS EQUIPMENT IN A RESIDENTIAL AREA IS LIKELY TO CAUSE INTERFERENCE IN WHICH CASE THE USER AT HIS OWN EXPENSE WILL BE REQUIRED TO TAKE WHATEVER MEASURES MAY BE REQUIRED TO CORRECT THE INTERFERENCE.

SECTION I

THEORY OF OPERATION

GENERAL OPERATING DESCRIPTION

The PHI Model 96A V/F Preamp detects analyzer collector current and converts the output to a proportional frequency or amplified analog signal.

INPUT AMPLIFIERS

The voltage-to-frequency preamp receives an input current from the collector of the analyzer electron multiplier. Referring to the block diagram, Figure 1, the first stage converts the current to a signal voltage with a bandwidth of about 10 kHz. An output connector is available at this point which is marked LOCK-IN. The signal is also further amplified by the second stage after which it goes to a voltage-to-frequency converter.

V/F CONVERTER

This stage produces a pulse train, whose frequency is linearly proportional to the input and whose maximum frequency is 1 MHz. The pulse train is sent to a one-shot multivibrator which reduces the width of each pulse from 500 nanoseconds to 100 nanoseconds. This is done so that when the pulses are sent through the pulse transformer at the 1 MHz rate, that this transformer has time to recover and consequently will not cause the secondary voltage to dc level shift.

ISOLATION TRANSFORMER

A pulse transformer is used as the coupling medium to shift the signal reference to ground level from the multiplier voltage level.

OUTPUT AMPLIFIERS

After the signal is transmitted through the transformer, it is sent through a second one-shot multivibrator which widens the pulses back to their original 500 nanoseconds width. This is done because the next stage frequency-to-voltage converter requires it for proper operation. The FREQ OUT connector is the output signal of the second one-shot multivibrator. The frequency-to-voltage converter reconverts the pulse train back to an analog signal. The last stage in the unit takes the voltage signal from the frequency-to-voltage converter and amplifies once more before it is sent to the ANALOG OUT connector.

POWER SUPPLIES

A floating power supply is used to power the stages before the pulse transformer. Primary power for this floating supply and the voltages to run the ground referenced circuitry are obtained from the PHI Model 32-100 Multiplier Supply.




Figure 1 Block diagram, voltage-to-frequency amplifier.

SECTION II

CALIBRATION AND MAINTENANCE

CALIBRATION SETUP

This calibration procedure provides for bench calibration of the V/F Preamp using standard test equipment.

/ WARNING /

Electric Shock Hazard

High voltages are present when the preamp is connected to the PHI Model 32-100 Multiplier Supply. Disconnect cables before calibrating. Refer calibration to qualified personnel.

Test Equipment Needed

The following test equipment is needed to calibrate the preamp:

- DVM (3 1/2 digit).
- Universal Power Supply.
- Test Cables.
- Oscilloscope.
- Frequency Counter.
- Universal Current Source.
- Special Test Cover (provides access holes for potentiometers and cables).

Visual Check

Visually inspect printed circuit board for open or shorted conductors.

Visually inspect all wiring inside the box for possible shorts or opens.

Check for tightness of mounting of all components.

ALIGNMENT, CALIBRATION AND TESTING

All measurements are ± 10 percent unless stated otherwise.

Power Supply Check

Connect the Universal Power Supply with ± 15 V and +5 V to the POWER-IN connector J1 on the preamp box. Turn the power switch ON. With the dvm, measure the voltage as follows:

IC AR1 - Pin 7 + $14\text{ V} \pm 1\text{ V}$

AR1 - Pin 4 - $14\text{ V} \pm 1\text{ V}$

Note: Common on AR1 - Pin 3
Short C9 with a clip lead.

AR1 Offset Adjustment

With the dvm, measure the voltage at TP2. Adjust R23 for $0\text{ V} \pm 1\text{ mV}$.

V/F Adjustment

1. With the oscilloscope, measure at FREQ OUT.
 - Adjust R18 for approximately 0 Hz.
 - Adjust R15 for 0.45 V , $\pm 1\text{ mV}$ on TP1.
 - Adjust R37 for 0 V , $\pm 1\text{ mV}$ on ANALOG OUTPUT.
2. Connect the Universal Current Source to the COLLECTOR connector J4 and the HV/POS IN connector J6. Set the current source for -500 nA .

3. With dvm, measure the voltage at TP2. It should be approximately 10 V. Install the test cover.
4. With the Frequency Counter, measure at FREQ OUT. Adjust R8 for 1 mHz, ± 100 Hz.
5. With scope, check for 500 nsec ± 100 nsec pulse at FREQ OUT.
6. Adjust R28 for 10 V, ± 50 mV on ANALOG OUT.
7. Set Current Source to -5 nA.
 - Adjust R18 for 10 kHz on FREQ OUT.
 - Adjust R37 for 0.1 V on ANALOG OUT.
8. Set Current Source to -500 nA. Repeat steps 4-7.
9. Disconnect the current source. Adjust R23 for 3 kHz, ± 10 Hz at FREQ OUT.
10. Remove test cover, remove clip lead and install cover.





**CABLE INTERCONNECT DATA
COMPUTER GENERATED**

[illegible]

PERKIN-ELMER		TITLE		SHEET: 1		1/27/88		PART NUMBER		ENG. DWG NO/REV		EFF DATE	
PHYSICAL ELEC. DIV.		BOARD ASSY-0 TO F						616800				1B	
ITEM	PART NO.	QTY	U/M	DESCRIPTION				DRAWING NO.	REFERENCE	FROM TO			
001	616799	1.000	EA	PC BOARD-0 TO F									
002	603716	1.000	EA	IC-DC/DC CONVERTER, ANALOG# 244					U4				
003	177228	2.000	EA	IC-OP AMP, JFET, INPUT, LF356N					AR2,3				
004	472038	1.000	EA	IC-OP AMP, 72301P					AR1				
005	616777	2.000	EA	IC-V/F/V CONVERTER, AD650KN					U1,3				
006	605912	4.000	EA	XSTR-NPN, 15V, .2A, 2N2369A					G1,2,3,4				
007	372075	1.000	EA	IC-TTL, MONOSTABLE MV, 74121					U2				
008	602791	2.000	EA	DIODE-SCH BARR, 15V, .02A, AS2811					CR7,8				
009	603955	7.000	EA	DIODE-HS SW, 75V, .25W, 1N914					CR2,3,4,4A,4B,5A,5B				
010	273016	1.000	EA	DIODE-ZENER, 10V, 400MW, 1N758A					CR5				
011	605911	1.000	EA	CORE-TOROIDAL, 768T188-3E2A									
012	373054	2.000	EA	CHOKE-4.7MH					L1,2				
013	602668	1.000	EA	CAP- 33UF, 10VDC, DP TANT					C19				
014	602136	10.000	EA	CAP- .22UF, 50V, ML CER RDL					CF'S, C13, 16, 17, 18				
015	602028	1.000	EA	CAP- .10UF, SKVDC, FILM TUB					C9				
016	373045	3.000	EA	CAP- .15PF, 1KV, CER DISC					C1, 12, 15				
017	171636	1.000	EA	CAP- .002MF, 6KV, CER-DISC					C3				
019	472043	1.000	EA	CAP- 1 UF, 50V, ML CER, RDL					C8				
020	616776	2.000	EA	CAP- 56 PF, 1X, 100V, ML CER					C6, 14B				
021	371015	3.000	EA	CAP- .001MF, 1KV, CER					C7, 11, 14				
022	477157	4.000	EA	CAP- .1 UF, 50V, ML CER RDL					C5, 8A, 10, 14A				
023	174085	1.000	EA	CAP- 5PF, 1KV, CER DISC					C2				
025	373057	1.000	EA	POT- 10 K, .52, 20T, TRMR-T/A					R23				
026	177402	2.000	EA	POT- 20K, .25W, 20 TURNS, TRIM T/A					R8, 15A				
027	607420	2.000	EA	POT- 2K OHM, .25W, 21T, 10PPM, T/A					R6, 27				
028	279285	1.000	EA	RES- 15 N OHM, .25W, 5%, RCO7					R12				
029	279205	4.000	EA	RES- 100 OHM, .25W, 5%, RCO7					R21, 32, 32A, 33				
030	279223	1.000	EA	RES- 2.2K OHM, .25W, 5%, RCO7					R15				
031	279229	2.000	EA	RES- 1.5K OHM, .25W, 5%, RCO7					R7, 17				
032	602147	1.000	EA	RES- 3.6K OHM, .25W, 5%, RCO7					R18A				
033	603597	1.000	EA	RES- 20K OHM, .12W, 1X, RN55D					R31				
034	603599	1.000	EA	RES- 1M OHM, .12W, 1X, RN55D					R20				
035	604965	2.000	EA	RES- 1.02K OHM, .12W, 1X, RN55D					R22, 24				
036	604347	1.000	EA	RES- 49.9K OHM, .12W, 1X, RN55D					R25				
037	603131	1.000	EA	RES- 100 OHM, .12W, 1X, RN55D					R26				
038	604969	1.000	EA	RES- 36.5K OHM, .12W, 1X, RN55D					R19				
039	609303	2.000	EA	RES- 249K OHM, .12W, 1X, RN55D					R14, 18B				
040	679263	1.000	EA	RES- 360K OHM, .25W, 5%, RCO7					R16				
041	604345	1.000	EA	RES- 5.74K OHM, .12W, 1X, RN55D					R13				
042	603936	1.000	EA	RES- 243 OHM, .12W, 1X, RN55D					R36				
043	604967	1.000	EA	RES- 9.76K OHM, .12W, 1X, RN55D					R35				

DATE NO.	605978	REV.	1	REV.	0
		REVISIONS			
ZONE	REV.	DESCRIPTION			
	A	RELINSE PER. TO 323			
		SCHEDULE 602013			
	B	SEE ECO 4470			
	C	HOW CHANGE			
		SEE ECO # 7545			



NOTES:

- ① MARK APPLICABLE INFORMATION ON ID. LABEL, ITEM 15, PER ENGRG SPEC. 60328 AND AFFIX APPROX AS SHOWN (COMP. ID. NO. ATG).
- ② WIRES ENDING AT THIS CONNECTOR ARE IN TWISTED PAIRS (ONE RED, ONE BLACK) FROM EACH BNC CONNECTOR.

SEE DETACHED PART LIST

SEE DETACHED PART LIST		THIS DRAWING CONTAINS AND IS THE PROPERTY OF PERKIN ELMER CORPORATION. IT IS TO BE USED FOR REPRODUCTION PURPOSES ON THE BASIS OF THE FOLLOWING AGREEMENT: THIS DRAWING IS TO BE USED FOR THE MANUFACTURE OF THE ITEM IDENTIFIED BY THE PART NUMBER AND TITLE HEREON. IT IS NOT TO BE REPRODUCED OR USED FOR ANY OTHER PURPOSE WITHOUT THE WRITTEN PERMISSION OF PERKIN ELMER CORPORATION.		TITLE MODEL 96A - V/F REF.FM		PERKIN ELMER CORPORATION Physical Electronics Division 3500 Central Expressway Santa Clara, California 95051		DRAWING NO. 60597 SHEET 1 OF 1	
DATE 10/25/71 REV 000 DESIGNED BY DRAWN BY CHECKED BY APPROVED BY		DATE 10/25/71 DATE 10/25/71 DATE 10/25/71 DATE 10/25/71		FINISH MATERIAL		SIZE D		SHEET 1 OF 1	

PERKIN-ELMER		TITLE SHEET: 1 12/16/87		PART NUMBER		ENG. DWG NO/REV		EFF DATE	
PHYSICAL ELEC. DIV.		MOD 96A-V/F PREAMP		605978		10		FROM TO	
ITEM	PART NO.	QTY	UM	DESCRIPTION	DRAWING NO.	REFERENCE			
0001	605443	1.000	EA	HOUSING-V/F PREAMP;MODEL 96A	TN U/616800				10188
0002	605249	1.000	EA	BOARD ASSY-V TO F					
0002	616800	1.000	EA	BOARD ASSY-V TO F					
0003	601542	3.000	EA	CONN-COAX,SHV,CHAS RCPT 1704-1		J4,5,6		10188	
0004	171289	3.000	EA	CONN-COAX,ENC,CHAS,RCPT,31-221		J2,3,7			
0005	601834	1.000	EA	CONN-CIRC,CHAS RCPT,5P,MS3102A		J1			
0006	541806	4.000	EA	SCREW-FNM, 4-40X .375,SST PHH					
0007	542104	4.000	EA	NUT-HEX,SMALL PATTERN,4-40 SST					
0008	512104	4.000	EA	WASHER-INT,TOOTH,#4 ,SST					
0009	602123	1.000	EA	STDF-INSUL,TUR,6-32 INTNL THD		TB1			
0010	541904	2.000	EA	SCREW-FNM, 6-32X .250,SST,PHH					
0011	512106	2.000	EA	WASHER-INT,TOOTH,#6 ,SST					
0012	471072	1.000	EA	SPCR-RND,.38 OD,6-32,.5 LG,CER		TB2			
0013	501000	4.000	EA	LUG-PLAIN SOLDER,3/8 STUD					
0014	546307	6.000	EA	SCREW-FLM, 6-32X .437,SST,PHH	PHI--192				
0015	479092	1.000	EA	LABEL-UNIT ID,1.38 X 2.75,POLY					
0016	501106	1.000	EA	LUG-LOCKING SOLDER,#6 STUD					
0017	541903	1.000	EA	SCREW-FNM, 6-32X .187,SST,PHH					
0018	606002	4.000	EA	SPCR-HEX,.25AF,6-32,.75LG,NYL					
0019	606001	4.000	EA	SCR-FLM,6-32X.375,NYL/MET CORE					
0020	604083	4.000	EA	SCR-RHD,6-32X.25,NYL/MET CORE					
0021	177059	1.000	EA	CONN-PLUG,CRIMP TYPE,2 SKT					
0022	177271	1.000	EA	CONN-PLUG,CRIMP TYPE,5 SKT					
0023	377092	1.000	EA	CONN-HSG,6 POSN,W/LKG RAMP					
0024	177070	13.000	EA	TERM-CRP CONN CONT,FOR S/S HSG					
0025	171636	1.000	EA	CAP-.002MF,6KV,CER-DISC		C1			
0026	279228	1.000	EA	RES- 4.7 K, .25W, 5% RC07		R1			
0027	279252	2.000	EA	RES- 100K OHM,.25W, 5% RC07		R2,3			
0028	502201	.900	FT	WIRE-STRD,22GA,600V,IPVC,BLK					
0029	502202	.300	FT	WIRE-STRD,22GA,600V,IPVC,BRN					
0030	502203	1.000	FT	WIRE-STRD,22GA,600V,IPVC,RED					
0031	502204	.300	FT	WIRE-STRD,22GA,600V,IPVC,ORN					
0032	502205	.300	FT	WIRE-STRD,22GA,600V,IPVC,YEL					
0033	502208	.300	FT	WIRE-STRD,22GA,600V,IPVC,VIO					
0034	602076	.700	FT	WIRE-STRD,20GA,9K0VDC,TFL,RED					
0035	605980	.000	EA	FUNC SCHEM-V TO F MODEL 96A	IN U/616795				
0036	605981	.000	EA	TEST AND CAL PROC-MODEL 96A				10188	10188
0036	616795	.000	EA	TEST & CAL PROC-MODEL 96A					
0037	605982	.000	EA	PROD SPEC-MODEL 96A V TO F					
				9036 END OF REPORT 9036					

DETACHED PARTS LIST

PERKIN-ELMER
Physical Electronics

TITLE: BOARD ASSEMBLY
V/F PREAMP

INACTIVE - FOR SERVICE
USE ONLY - SUPERSEDED

PART NO. 605349
REV F

ITEM NO.	PART NUMBER	QUANTITY PER ASSEMBLY	DESCRIPTION	7 REF. DESIGNATIONS / NOTES
1	605348	1.	PC BD - V/F/V PREAMP	U4
2	603716	1.	IC-DC/DC CONVERTOR ANALOG 944	AR2, 3
3	177228	2.	IC-OP AMP/FET INPUT LF 356N	ARI
4	472038	1.	IC-OP AMP LM 301P	U1, 3
5	604970	2.	IC-V/F/V CONVERTOR INTECH A8404	Q1-4
6	605912	4.	TRANSISTOR - 2N2369A	U2
7	372075	1.	IC-TTL MONOSTABLE MV 74121	CR7, 8
8	602791	2.	DIODE-SCH BAR, 15V, 0.2A, AS2811	CR2-4
9	603955	3.	DIODE-H, SP, SW, 75V, 25W IN914	CR5
10	273016	1.	DIODE-ZENER, 10V, 400MW IN758A	
11	605911	1.	CORE-TOROIDAL, 768788-3E2A	L1, 2
12	373054	2.	CHOKE 4.7MH	C19
13	602668	1.	CAP-33MFD, 10VDC, DPD TANT.	C13, 16-18, C.F.(6)
14	602136	10.	A -.22MFD, 50VDC, MONO CER	C9
15	602038	1.	-.10MFD, 5000VDC, FILM TUBE	C1, 12, 15
16	373045	3.	-.15PFD, 1KV, CER. DISC.	C3
17	171636	1.	-.002MFD, 6KVDG, CER. DISC	C7
18	273021	1.	-.330PFD, 1KV, CER. DISC	C8, 10
19	472043	2.	-.1UF, 50V, CER. RDL	C6, 14
20	177256	2.	-.33PFD, 300N, DPD MICA	C11
21	171455	1.	-.005MFD, 100V, CER. DISC	C5
22	273002	1.	↓ -.01MFD, 50V, CER. RDL	C2
23	174085	1.	CAP-.5PFD, 1KV, CER. DISC	R18
24	472033	1.	POT-1K, .25W, 20T, TRIM, T/A	R23
25	373057	1.	POT-10K, .5W, 20T, TRIM, T/A	

PREPARED: *[Signature]* DATE: 12-10-71
 CHECKED: *[Signature]* DATE: 12-15-71
 APPROVED: *[Signature]* DATE: 12-15-71
 ADD: ☐

PERKIN-ELMER
Physical Electronics


PART NO. 605349 REV E

SHEET 3 OF 3


[illegible]


REVISIONS				
ZONE	REV.	DESCRIPTION	DATE	APPROVED
	A	RELEASE PER ECO 3123 SUPERSEDES 602013	12-5-81	AT
	B	SEE ECO 3256 WAA	4-22-82	AT
	C	HDW CHG	3-15-83	AT
	D	SEE ECO 5174 KWA	8-13-83	AT
	E	HDW CHG	8-20-85	AT
	F	SEE ECO # 7545	1-25-86	TLL

00761600 FBI FPO 7545
USE CIVL - SUPERSEED
EQUIV - POL SHOOT



DETAIL A (T1)
SCALE: NONE



ASSY 6053494 REV 
V TO F MADE IN USA

NOTES:

1. FABRICATE & ASSEMBLE PER PERKIN-ELMER PHI, SPECIFICATIONS.
2. MARK APPLICABLE REVISION LEVEL IN SPACE PROVIDED.
3. COVER WIRE WITH TUBING AND SHRINK COMPLETELY BEFORE WINDING ON CORE.
4. CHOKES L1 & L2 MUST BE INSTALLED WITH THEIR WINDINGS NOT CONTACTING GROUND PLANE ON THE PC BOARD.

SEE DETACHED PART LIST

[illegible]