


Before



After



### Description

With the exchange of your existing 04-725 (or 04-727, 07-763) sample intro arm, RBD can provide you with a magnetic load-lock system at a fraction of the cost of other magnetic load-lock upgrades.

### Advantages

- Better vacuum during sample loading
- System won't vent if the probe is bumped
- No Teflon seals and spacers to replace
- No expensive intro switch assembly to break
- 2.75" CF Mounting
- Magnetic switch compatible with Perkin-Elmer Auto Valve Controller
- Compatible with standard sample mounts