



RBD Instruments, Inc.
2437 NE Twin Knolls Dr.
Suite 2`
Bend, Oregon, USA

+1 541-330-0723

<https://rbdinstruments.com/products/vb.html>

VB FAQ

Is the length of the IRB-600 Emitter customizable?

Yes. The “standard” maximum length is 6" (4" beyond the standard length). Customization requests for an extended IRB-600 length that is longer than 6" will be evaluated on a case-by-case basis with respect to feasibility.

What is the average lifetime of the IRB-600 emitter?

The average lifetime is 5,000 hours. This can be affected by how hot the emitter is run and how close it is to other surfaces (in terms of reflected energy).

What is the standard length of the cable between the BC-3 and the IRB?

The standard cable length is 15 feet.

What is the maximum length of the cable between the BC-3 and the IRB?

The maximum cable length is 15 feet. This is the standard length.