



microCMA Dimensions Form

Name: _____

Company: _____

Date: _____

Comments: _____

Signature: _____

By signing this I agree that the dimensions entered on this form meet my specifications.

Flange to target distance _____ mm

+ Z-translator fully compressed length

| | | |
|---------------------------------------|---------------------------------------|--|
| 66 mm <i>(50 mm stroke)</i> | 70 mm <i>(75 mm stroke)</i> | 74 mm <i>(100 mm stroke)</i> |
| Custom _____ mm | | |

- End of optics to target 4 mm

= microCMA length _____ mm *(standard length 283 mm)*

Vacuum Chamber

